



*Your supplier of high-pressure laboratory instruments  
and advanced geotechnical testing equipment*

# HYDRATEST

## Gas Hydrate Triaxial Compression Test System

*Triaxial testing of sediments containing methane hydrates under extreme low-temperature  
and high-pressure conditions.*



## What is HYDRATEST?

- Servo-controlled triaxial testing system for gas hydrate-bearing sediments
- Operates under high pressure and sub-zero (below freezing) temperatures
- Investigates mechanical properties of Methane Hydrate-Bearing Sediments (MHBS)
- Simulates in-situ conditions: temperature, pore & confining pressure, hydrate saturation
- Full hydrate formation and dissociation cycles
- Accurate measurement of Young's modulus, cohesion, internal friction angle
- Automated control and data acquisition via Floxlab proprietary software

## KEY SPECIFICATIONS

**100 kN**

Axial Load Capacity

**70 MPa**


Max Confining & Pore Pressure

**-10 → +60 °C**

Operating Temperature

**50 mm**

Specimen Diameter



## TEST PRINCIPLE

- 1 Sand specimen initially saturated with 10% water, sealed in a Teflon sleeve
- 2 Placed in triaxial cell; water injected by syringe pump until target saturation
- 3 Methane injected at controlled flow rate under confining pressure
- 4 Constant pore pressure for 24 h → CH<sub>4</sub> dissolution in pore water
- 5 Temperature lowered into hydrate-stability range → MH formation
- 6 Axial load applied; strain, force, pore pressure recorded
- 7 Temperature raised → MH dissociation; released CH<sub>4</sub> volume measured

## TESTS PERFORMED

- Triaxial tests on gas hydrate sediment specimens
- Hydrate formation tests
- Hydrate dissociation tests
- Drained & undrained compression tests
- Creep & strength characterization of MHBS
- Controlled cooling, heating & thermal cycling

## MEASURED RESULTS

- Young's modulus
- Poisson's ratio
- Mohr circle envelope
- Axial stress & strain
- Diametral strain
- Pore pressure
- Confining pressure
- Specimen temperature
- Water pressure & volume
- Produced CH<sub>4</sub> gas volume



Complete HYDRATEST system

Deep-sea conditions

Cryogenic capability

Full automation

## SCOPE OF SUPPLY

- Servo-controlled hydraulic load frame (100 kN)
- High-pressure triaxial loading cell
- In-vessel axial & diametral deformation sensors
- In-vessel strain-gauge load cell
- Thermostatic bath with cooling / heating jacket
- Confining pressure intensifier
- Water syringe pump
- Methane gas syringe pump
- Produced gas measurement unit (BPR + flowmeter)
- Gas booster (up to 700 bar)
- Computer workstation with proprietary Floxlab software

## LOADING SYSTEM

<b>Servo-Hydraulic Load Frame</b>	100 kN • stiffness 70 kN/mm • piston stroke 50 mm
<b>Triaxial Cell</b>	100 kN axial • 70 MPa • -10 to +60 °C • ø50 mm specimen

## PRESSURE & FLUID CONTROL

<b>Confining Pressure Intensifier</b>	70 MPa • accuracy 0.25% FS • 250 cc • 20-litre reservoir
<b>Water Syringe Pump</b>	70 MPa • 500 cc • flow 0.0001–100 cc/min
<b>Methane Gas Syringe Pump</b>	70 MPa • 500 cc • flow 0.0001–100 cc/min

## IN-VESSEL INSTRUMENTATION

<b>Diametral Extensometer</b>	range 2.5 mm • linearity 0.3% FS • -10 to +60 °C
<b>Axial LVDT Transducer</b>	±5 mm • linearity 0.15% FS • 3 LVDTs at 120°
<b>Strain-Gauged Load Cell</b>	100 kN • accuracy ±0.5% • sensitivity ≈ 2 mV/V

## THERMAL & GAS MANAGEMENT

<b>Thermostatic Bath &amp; Cooling Jacket</b>	-10 to +60 °C • resolution 0.1 °C • 10 L/min • 4.5 L
<b>Produced Gas Measurement Unit</b>	BPR up to 10,000 psi • 0–1,000 sccm • stainless steel
<b>Gas Booster (GB 700)</b>	inlet 15–200 bar • outlet up to 700 bar • ratio 1:75

## CONTROL & ACQUISITION

<b>FloXlab Proprietary Software</b>	real-time monitoring • automated sequences • synchronized logging • standardized reports
-------------------------------------	--



# GET IN TOUCH

*We look forward to discussing your project.*

■ **Address**

23 rue du Port • 92000 Nanterre — France

● **Phone**

+33 (0)1 81 93 12 85

■ **Email**

[contact@floxlab.com](mailto:contact@floxlab.com)



SCAN ME

[www.floxlab.com](http://www.floxlab.com)