

FLOXLAB AUTO-VALVE FEATURES



Model	Pressure (psi)	Temperature (°C)	Wetted parts	Seal	Connection
TWO-WAY ON/OFF VALVES					
SS03-AT	3,000	80	Stainless steel	Viton	1/4"LOK
SS03-HT	3,000	160	Stainless steel	Viton	1/4"LOK
SH03-AT	3,000	80	Hastelloy	Viton	1/4"LOK
SH03-HT	3,000	160	Hastelloy	Viton	1/4"LOK
SS10HF-AT	10,000	80	Stainless steel	Viton	1/4"SLP
SS10HF-HT	10,000	160	Stainless steel	Viton	1/4"SLP
SH10HF-AT	10,000	80	Hastelloy	Viton	1/4"SLP
SH10HF-HT	10,000	160	Hastelloy	Viton	1/4"SLP
SS10-AT	10,000	80	Stainless steel	Viton	1/8"LP
SS10-HT	10,000	160	Stainless steel	Viton	1/8"LP
SH10-AT	10,000	80	Hastelloy	Viton	1/8"LP
SH10-HT	10,000	160	Hastelloy	Viton	1/8"LP
SS10VF-AT	10,000	80	Stainless steel	Viton	3/8"SLP
SS10FV-HT	10,000	160	Stainless steel	Viton	3/8"SLP
SS20-AT	20,000	80	Stainless steel	Viton	1/4"HP
SS20-HT	20,000	160	Stainless steel	Viton	1/4"HP
SH20-AT	20,000	80	Hastelloy	Viton	1/4"HP
SH20-HT	20,000	160	Hastelloy	Viton	1/4"HP
SS30-AT	30,000	40	Stainless steel	Polyurethane	1/4"HP
SH30-AT	30,000	40	Hastelloy	Polyurethane	1/4"HP
THREE-WAY 4-POSITION VALVES					
DS10HF-AT	10,000	80	Stainless steel	Viton	1/4"SLP
DS10HF-HT	10,000	160	Stainless steel	Viton	1/4"SLP
DH10HF-AT	10,000	80	Hastelloy	Viton	1/4"SLP
DH10HF-HT	10,000	160	Hastelloy	Viton	1/4"SLP
DS10-AT	10,000	80	Stainless steel	Viton	1/8"LP
DS10-HT	10,000	160	Stainless steel	Viton	1/8"LP
DH10-AT	10,000	80	Hastelloy	Viton	1/8"LP
DH10-HT	10,000	160	Hastelloy	Viton	1/8"LP
DS10VF-AT	10,000	80	Stainless steel	Viton	3/8"SLP
DS10FV-HT	10,000	160	Stainless steel	Viton	3/8"SLP
DS20-AT	20,000	80	Stainless steel	Viton	1/4"HP
DS20-HT	20,000	160	Stainless steel	Viton	1/4"HP
DH20-AT	20,000	80	Hastelloy	Viton	1/4"HP
DH20-HT	20,000	160	Hastelloy	Viton	1/4"HP
DS30-AT	30,000	40	Stainless steel	Polyurethane	1/4"HP
DH30-AT	30,000	40	Hastelloy	Polyurethane	1/4"HP

Note: Requires 80-115 psi air